



SELDI ANALYSIS

The Need for the Review and Understanding of SELDI/MALDI Data Prior to Analysis (Analyzer Beware)

William E. Grizzle, O. John Semmes, William Bigbee, Liu Zhu,
Gunjan Malik, Denise K. Oelschlager, Barkha Manne, Upender
Manne, University of Alabama at Birmingham, Eastern Virginia
Medical School, University of Pittsburgh Medical Center



SELDI ANALYSIS

WINDOW INTO DISEASES

Plasma Serum Urine

Effusions

Mucous Saliva Fecal Matter

Tissue (PAP; Urinary Sediment)

Bile CSF Sweat

SELDI ANALYSIS







Illustration of SELDI Time-Of-Flight (TOF) Mass Spectrometry.
(Modified with permission from Ciphergen Biosystems, Inc.)





**Contents of Duct
Including
Products of
Dying Cells and
Living Cells e.g.,
PSA**




**Products of Dying
Cancer Cells Collect in
Interstitial Space and are
Absorbed into Vascular
and Lymphatic Vessels**



**Products from Dying
Cells and Duct
Contents e.g., PSA**






ONCOFETAL TUMOR ANTIGENS AND METABOLITES



CEA
TAG-72
CA125
MUC-2
LEWIS Y

MOLECULES AND
METABOLITES OF NORMAL
TISSUES - ACTIN







SELDI ANALYSIS

JUNK



**Statistician
Bioinformatician**



JUNK



SELDI ANALYSIS

PROBLEMATIC ISSUES IN ANALYSIS

Experimental Design

Patient

Sample

Protein Chip

Spectra

Analytical Approach



SELDI ANALYSIS

Experimental Design

Selecting Cases and Controls

Collecting and Processing Samples

Performing Assays Without Bias

Selecting Optimal Approach to
Analysis

Avoiding Over-Analysis



SELDI ANALYSIS

PATIENT

Groups Comparable

Sites

Racial/Ethnic Balance

Homeostatic Balance

No Bias



SELDI ANALYSIS

PATIENT

Usually Comparing Disease vs. Control

or

Disease A vs. Condition B vs. Normal



SELDI ANALYSIS

PATIENT

Control Definition

Is Disease Absent?

Is There Bias In The Controls?



SELDI ANALYSIS

SAMPLE

Type

Collection

Processing

Storage

Transfer




SELDI ANALYSIS

SAMPLE


Type

Serum versus Plasma

Sensitivity- Can Products of
Tumors Be Detected



LNCAP CELL LYSATES USING WCX2 ARRAYS



100 Cells per
ml

100 Cells per
ml

1000 Cells per
ml

1000 Cells per
ml

10,000 Cells per
ml

10,000 Cells per
ml

10⁵ Cells per
ml

10⁵ Cells per
ml



SELDI ANALYSIS

SAMPLE

Collection

How?

Stress

Container-e.g., multiple
anticoagulants; thus
all Plasmas are not
the Same



SELDI ANALYSIS

SAMPLE

Processing

Time from Collection to
Freezing (Too Restrictive?)

For Consistent Results, Robotic
Processing Is Required



SELDI ANALYSIS

SAMPLE

Processing

Removal of Proteins Present in
Large Concentrations
(e.g., Albumin) May Also
Remove Peptides Being
Carried by Removed
Proteins



SELDI ANALYSIS


SAMPLE

Storage

Length

Temperature

QC Sample Spectra



Aliquot of Sample (#20) Stored at -80C

Aliquot of Sample (#20) Transferred from -80C and Stored for 3 Months at -20C

Aliquot of Sample (#20) Stored at -20C for 5 Months

Aliquot of Sample (#20) Stored at -20C for 7 Months

Aliquot of Sample (#20) Stored at -20C for 8 Months

Second Aliquot of Sample (#20) Stored at -20C for 8 Months

Second Aliquot of Sample (#20) Stored ONLY at -80C for Ten Additional Months




SELDI ANALYSIS

Transfer

Freeze-Thaw Cycles

Quantity (Triplicates and Repeat-300
mcl)

Decrease in Protein Intensity in Human Serum due to repetitive Freeze-Thaws



Decrease in Protein Intensity in Human Serum (#6) due to repetitive Freeze-Thaws




Figure1: The graph above shows the decrease in protein intensities in three different proteins in one human serum sample. The legend indicates the molecular weight of each protein. The specific intensity values are displayed above each of the bars.



SELDI ANALYSIS

- **Chip Type**






TABLE 1:

Old Designation	Current Chip	Biochemical Action of Surface Chemistry
IMAC3	IMAC30 (with hydrophobic barrier)	Bivalent metals can be attached to the chip. Proteins that bind to these divalent metals (e.g., Cu^{+2}) are bound by the chip.
WCX2	Same (CM10 mimics WCX2 but does not replace it)	This is a weak cation exchange chip. It contains negatively charged (anionic) carboxylate groups that will bind proteins with positively charged areas containing high numbers of lysine, arginine, and/or histidine amino acids.
H4	Same (C16 contains 16 CH_3)	The chip contains multiple chains of 16 methylene groups. This binds molecules that are hydrophobic.
SAX2	Q10 (with hydrophobic barrier)	Strong anion exchanger which is composed of quaternary ammonium groups that are charged positively. This chip will bind proteins/peptides with regions rich in acidic groups, especially regions of peptides high in aspartic and/or glutamic amino acids.
NP1 and NP2	NP20	General protein binding surface with binding of hydrophilic proteins.
PS1 and PS2	PS10 / PS20	Chip designed to bind capture molecules of choices e.g., antibodies, receptors, and nuclei acid binding proteins PS-1 (carbonyl diimidazole groups), PS-2 (epoxy groups). Also the PS-2 has a hydrophobic coating.
SENDID		Incorporates EAM into chip.

QC Sample Spectra at Different Biochips





SELDI ANALYSIS

Protein Chip

- High Concentration Proteins May Block Binding of Low Concentration Proteins:
10,000 Ci of 5500 D Protein vs. 10 Ci of 7500 D Protein both with Same Binding Characteristics




SELDI ANALYSIS

Protein Chip

- All Proteins/Peptides Bound to Chip May Not Be Released/Ionized.

QC Sample Spectra at Multi-times Reading



FIRST READ A

FIRST READ B

FIRST READ C

4TH READ B


10TH READ B

14TH READ B

SELDI ANALYSIS

Spectrum

- “Directed” and “Non-Directed” Approaches to Begin Spectral Analysis
- Directed= Peak at 5500 Same As Peak at 5507 Based on Resolution $+0.2\%$





SELDI ANALYSIS

Spectrum


- Primary Peaks- Disease Has Unique or Larger Peak than Non-Disease; Thus, the Disease Produces a Molecular Product.
- Secondary Peak-Disease Causes a Decrease in Molecular Species Normally Present via Change in Metabolism or Excretion and/or Shutdown in Production



SELDI ANALYSIS

Spectrum

- Components of Spectra at Molecular Weights of Less Than 20,000 May Represent Metabolites of Proteins/Peptides Rather Than Intact Proteins/Peptides






SELDI ANALYSIS

Spectrum

- Peaks May Not Provide Independent Information: For Example the Peak at 5500 D May Be A Metabolite of the Peak at 7500

7500 → 5500






SELDI ANALYSIS

Spectrum

- A High Concentration Protein May Prevent Identification of Low Concentration Protein: 1200 Ci of 5500 D Protein vs. 100 Ci of 5510 D Protein Even with Different Binding Characteristics to Same Chip






SELDI ANALYSIS

Spectrum

- All Areas of the Spectrum Are Not The Same
 - Molecular Weights of Less Than 2000
 - No Standards; Noise; Contamination
 - Weights of Greater Than 50,000
 - Proteins of High Concentration






SELDI ANALYSIS

Spectrum

- How Variable Is The Peak Location and/or Amplitude When the Same Sample Is Run On the Same Day on the Same Machine?
On the Next Day on the Same Machine?
On a Different Machine?






SELDI ANALYSIS

Spectrum

- Eastern Virginia Medical School; UAB; U of Texas San Antonio; U of Pittsburgh Medical Center; Johns Hopkins Medical Center; Uniformed Health Services
- ***All Were Able To Standardize Their Machines and To Obtain Comparable Data on 14 Cancer and 14 Non-Cancer Cases***



SELDI ANALYSIS
